

Aminomethyltransferase Polyclonal Antibody

Catalog number: 10633-1-AP

Size: 20 µg/150 µl

Source: Rabbit

Isotype: IgG

Synonyms:

AMT; Aminomethyltransferase,

AMT, GCE, GCST, GCVT, NKH

Applications

Tested applications: ELISA
Species specificity: Human, Mouse, Rat; other species not tested.
Calculated 44kd
Aminomethyltransferase
MW:
Observed
Aminomethyltransferase
MW:

Application key: WB = Western blotting, IHC = Immunohistochemistry, IF = Immunofluorescence, IP = Immunoprecipitation

Immunogen information

Immunogen: Ag0990
GenBank accession number: BC007546
Gene ID (NCBI): 275
Full name: Aminomethyltransferase

Product information

Purification method: Antigen affinity purification
Storage: PBS with 0.1% sodium azide and 50% glycerol pH 7.3. Store at -20°C.