

Specification Sheet

PRODUCT: Water, USP Purified [Sterile Filtered, WFI Quality]

CATALOG NO.:	462918	LOT NO.:	TBD
MANUFACTURE DATE	TBD	EXPIRATION DATE:	12 Months from the date of manufacture
STORAGE CONDITIONS:	4 - 40°C	CONTAINER SIZE:	20L

Test Results

Test	Specification	Test Method	Result
Conductivity	Meets Requirements ≤ 1.3 µS/cm @ 25°C	USP <645> Bulk Water	TBD
Total Organic Carbon	≤ 500 ppb	USP <643> Bulk Water	TBD
Sterility	No Growth	USP <71>	TBD
Endotoxin	≤ 0.25 EU/mL	USP <85>	TBD

Not for direct human, animal or therapeutic use. For Research and Manufacturing Uses.

Manufacturing Site

Lifecycle Biotechnologies, 5858 Wright Drive, Loveland, Colorado 80538