

LDB2 Polyclonal Antibody

Catalog number: 11873-1-AP

Size: 45 µg/150 µl

Source: Rabbit

Isotype: IgG

Synonyms:

LDB2; CLIM 1, CLIM1, LDB 2,

LDB1, LDB2, LIM domain

binding 2, LIM domain binding

protein 2

Applications

Tested applications: ELISA
Species specificity: Human, Mouse, Rat; other species not tested.
Calculated LDB2 MW: 42 kDa
Observed LDB2 MW:

Application key: WB = Western blotting, IHC = Immunohistochemistry, IF = Immunofluorescence, IP = Immunoprecipitation

Immunogen information

Immunogen: Ag2460
GenBank accession number: BC034019
Gene ID (NCBI): 9079
Full name: LIM domain binding 2

Product information

Purification method: Antigen affinity purification
Storage: PBS with 0.1% sodium azide and 50% glycerol pH 7.3. Store at -20°C.