

## DCTN5 Polyclonal Antibody

Catalog number: 10182-1-AP


Size: 25 µg/150 µl

Source: Rabbit

Isotype: IgG

Synonyms:

DCTN5; DCTN5, Dynactin  
subunit 5, Dynactin subunit  
p25, p25, dynactin 5 (p25)



Immunohistochemical of paraffin-embedded human breast cancer using 10182-1-AP (DCTN5 antibody) at dilution of 1:50 (under 10x lens)

### Applications

Tested applications:	ELISA, IHC
Cited applications:	WB
Species specificity:	Human, mouse, rat; other species not tested.
Cited species:	Monkey
Calculated DCTN5 MW:	20 kDa
Observed DCTN5 MW:	
Positive IHC detected in	Human breast cancer
Recommended dilution:	IHC: 1:20-1:200

Application key: WB = Western blotting, IHC = Immunohistochemistry, IF = Immunofluorescence, IP = Immunoprecipitation

### Immunogen information

Immunogen:	Ag0234
GenBank accession number:	BC004191
Gene ID (NCBI):	84516
Full name:	Dynactin 5 (p25)

### Product information

Purification method:	Antigen affinity purification
Storage:	PBS with 0.1% sodium azide and 50% glycerol pH 7.3. Store at -20°C.