

Puritan dry transport systems.

GET ALL OF YOUR SPECIMEN COLLECTION DEVICES IN ONE PLACE! Puritan has added dry transport systems – made in the USA – to our specimen collection and transport device lines. There’s a wide selection of standard material and patented HydraFlock® and PurFlock Ultra® tips, and many handle types (some with molded breakpoints). All models come sterile or non-sterile. Custom orders and configurations are always welcome.



Safe, sterile specimen transport. No matter which swab(s) you choose, all of our transport systems include a round-bottom tube made from virgin polypropylene resin, which will securely hold the swab, maintaining a safe and sterile environment.

DRY TRANSPORT SYSTEMS WITH A VARIETY OF STANDARD TIPPED APPLICATORS

These are designed to assure integrity of the sample collection and transport. We work closely with our customers to provide a range of tip and handle configurations, with the tube, to fit specific situations and environments. Custom orders welcome.

Item No.	Item Description	Cap Color	Handle length x diameter	Tip length x diameter	Tube length x diameter	Packaging pack, weight
sterile						
25-806 1PA BT	Regular calcium alginate tip, polystyrene handle.	Blue	149.225 x 2.54 mm 5.875 x 0.1"	17.45 x 5.156 mm 0.687 x 0.203"	154 x 12 mm 6.062 x 0.472"	5/100/1,* 7 lbs.
25-806 1WA BT	Regular calcium alginate tip, wood handle.	Blue	147.625 x 2.54 mm 5.812 x 0.1"	16.281 x 5.563 mm 0.641 x 0.219"		
25-806 1PC BT	Regular cotton tip, polystyrene handle.	Pink	147.225 x 2.54 mm 5.875 x 0.1"	17.45 x 5.156 mm 0.687 x 0.203"		
25-806 1WC BT	Regular cotton tip, wood handle.	Pink	147.625 x 2.54 mm 5.812 x 0.1"	15.875 x 4.775 mm 0.625 x 0.188"		
25-806 1PD BT	Regular polyester tip, polystyrene handle.	Green	149.225 x 2.54 mm 5.875 x 0.1"	17.45 x 5.156 mm 0.687 x 0.203"		
25-806 1PR BT	Regular rayon tip, polystyrene handle.	White	149.225 x 2.54 mm 5.875 x 0.1"	17.450 x 5.156 mm 0.687 x 0.203"		
25-806 1WR BT	Regular rayon tip, wood handle.	White	147.625 x 2.54 mm 5.812 x 0.1"	15.875 x 4.775 mm 0.625 x 0.188"		
25-1406 1PF BT	Small foam tip, polystyrene handle.	Blue	152.298 x 2.489 mm tapered to 1.321 mm 5.969 x 0.098" tapered to 0.049"	13.894 x 3.581 mm 0.547 x 0.141"		
25-1506 1PF BT	Regular foam tip, polystyrene handle.	Blue	149.225 x 2.54 mm 5.875 x 0.1"	16 x 4.775 mm 0.63 x 0.188"		



Color coded for confidence, and convenience.
 The injection-molded cap gives you a precision fit for multiple handle styles. Printed product numbers and colored caps assist in inventory maintenance and correct configuration in the field and in the lab. All Puritan dry transport systems are made in the USA. Custom sizes, materials, and configurations are always welcome.

DRY TRANSPORT SYSTEMS WITH PATENTED HYDRAFLOCK® & PURFLOCK ULTRA® FLOCKED SWABS

Our patented HydraFlock® (product numbers coded with an H) and PurFlock Ultra® (code U) flocked swabs are available in different tip shapes and handle configurations. It's the processing of the tip fibers themselves that creates the characteristic rapid absorption and quick release of specimens.

HydraFlock Item No.	PurFlock Ultra Item No.	Item Description	Cap Color	Handle length x diameter	Tip length x diameter	Tube length x diameter	Packaging pack, weight
sterile	sterile						
25-3306-H BT	25-3306-U BT	Standard tip, polystyrene handle with breakpoint at 80mm.	Red	152.654 x 2.54 mm tapered to 1.676 mm 6.01 x 0.1" tapered to 0.066"	17.45 x 5.08 mm 0.687 x 0.2"	154 x 12 mm 6.062 x 0.472"	5/100/1, 7 lbs.
25-3316-H BT	25-3316-U BT	Mini-tip, polystyrene handle, no breakpoint.	Black	149.6 x 2.489 mm tapered to 1.245 mm 5.89 x 0.098" tapered to 0.049"	18.415 x 3.429 mm 0.725 x 0.135"		
25-3316-H EMB 80 BT	25-3316-U EMB 80 BT	Mini-tip, polystyrene handle with breakpoint at 80mm.	Green	149.606 x 2.489 mm tapered to 1.245 mm 5.89 x 0.098" tapered to 0.049"	18.415 x 3.429 mm 0.725 x 0.135"		
25-3316-H EMB 100 BT	25-3316-U EMB 100 BT	Mini-tip, polystyrene handle with breakpoint at 100mm.	Yellow	149.606 x 2.489 mm tapered to 1.245 mm 5.89 x 0.098" tapered to 0.049"	18.415 x 3.429 mm 0.725 x 0.135"		
25-3317-H BT	25-3317-U BT	Ultrafine tip, polystyrene handle with breakpoint at 100mm.	Orange	152.146 x 2.54 mm tapered to 0.914 mm 5.99 x 0.1" tapered to 0.036"	16.764 x 3.556 0.66 x 0.14"		
25-3318-H BT	25-3318-U BT	Micro ultrafine tip, polystyrene handle with breakpoint at 100mm.	White	152.146 x 2.54 mm tapered to 0.914 mm 5.99 x 0.1" tapered to 0.036"	7.925 x 3.556 mm 0.312 x 0.14"		
25-3406-H BT	25-3406-U BT	Large tip, polystyrene handle with breakpoint at 80mm.	Red	149.606 x 2.54 mm tapered to 2.286 mm 5.89 x 0.1" tapered to 0.09"	16.256 x 5.563 mm 0.64 x 0.219"		
25-3606-H BT	25-3606-U BT	Elongated tip, polystyrene handle with breakpoint at 30mm.	White	151.384 x 2.54 mm tapered to 1.702 mm 5.96 x 0.1" tapered to 0.067"	17.399 x 5.63 mm 0.685 x 0.219"		
25-3706-H BT	25-3706-U BT	Elongated tip, polystyrene handle with breakpoint at 100mm.	White	151.384 x 2.54 mm tapered to 1.702 mm 5.96 x 0.1" tapered to 0.067"	17.399 x 5.63 mm 0.685 x 0.219"		
25-3806-H BT	25-3806-U BT	Elongated tip, polystyrene handle, no breakpoint.	White	151.384 x 2.54 mm tapered to 1.702 mm 5.96 x 0.1" tapered to 0.067"	17.399 x 5.63 mm 0.685 x 0.219"		

*5 boxes per case, 100 tubes per box, 1 swab per tube



207-876-3311 • puritandiagnostics.com
 sales@puritanmedproducts.com
 Puritan Diagnostics LLC, An affiliate of Puritan Medical Products Company
 31 School Street, Guilford, Maine 04443-0149 USA
 ISO 9001:2008 ISO 13485:2003 CE



facebook.com/PuritanUSA
 twitter.com/PuritanUSA
 linkedin.com/company/PuritanUSA
 youtube.com/PuritanUSA