

CACNB1 Polyclonal Antibody

Catalog number: 13039-1-AP


Size: 20 µg/150 µl

Source: Rabbit

Isotype: IgG

Synonyms:

CACNB1; CACNB1, CACNLB1,
CCHLB1, CAB1



mouse brain tissue were subjected to SDS PAGE followed by western blot with 13039-1-AP(CACNB1 antibody) at dilution of 1:500

Applications

Tested applications:	ELISA, WB
Species specificity:	Human, mouse, rat; other species not tested.
Calculated CACNB1 MW:	598aa,66 kDa
Observed CACNB1 MW:	75kd
Positive WB detected in	Mouse brain tissue
Recommended dilution:	WB: 1:200-1:2000

Application key: WB = Western blotting, IHC = Immunohistochemistry, IF = Immunofluorescence, IP = Immunoprecipitation

Immunogen information

Immunogen:	Ag3693
GenBank accession number:	BC037311
Gene ID (NCBI):	782
Full name:	Calcium channel, voltage-dependent, beta 1 subunit

Product information

Purification method:	Antigen affinity purification
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3. Store at -20°C.