




Specifications

Electrical requirements	110VAC 60Hz single phase, 1700 watts
Chamber dimensions	17 3/4 x 11 1/2 x 18" (WxDxH)
Unit dimensions	24 x 22 x 30.5"(WxDxH)
Shipping dimensions	27 x 30 x 27" (WxDxH)
Temperature range	Ambient temperature to 300°C
Temperature uniformity	+/- 2°F
Temperature controller precision	0.1°F
Timer range	1 to 9999 minutes
Racks included	2 (additional racks available at extra cost)
Maximum rack capacity	7 (height between each rack is 2 inches)
Exhaust port	2" diameter x 1, adjustable
Shipping weight	190 Lb
Access port	2" diameter with cover, one
Certification	CE, UL/CSA (optional)

Standard package

Part description	Quantity	Part image
Ai 2.5 cu ft forced air oven	1 pc	
Stainless steel rack and brackets	2 sets	
110 volt 3-prong power cord	1 pc	
User's manual	1 pc	