



DETAIL B  
SCALE 2 : 1



Mass properties of water test

Density = 1.00000000 grams per cubic centimeter

Mass = 1.15887005 grams

Volume = 1.15887005 cubic centimeters

This Mold Is Owned by KDeb Manufacturing Corporation

DO NOT SCALE DRAWING		NAME		DATE		<b>BATESVILLE TOOLING &amp; DESIGN, INC.</b>	
TOLERANCES UNLESS OTHERWISE SPECIFIED: DIMENSIONS ARE IN INCHES		DRAWN BY		gsb		210 TOWER ROAD BATESVILLE, MS 38606	
TOLERANCES:		CHECKED BY				W.M. HARMON IND COMPLEX	
FRACTIONAL ± 1/64		APPROVED BY				(662)563-0235 FX (662)563-1663 PH	
ANGULAR: MACH ± 1/2 DEG		MATERIAL:		Polypropylene Homopolymer Ineos H12Z-00 or Equivalent		www.btd-inc.net	
TWO PLACE DECIMAL ±.005							
THREE PLACE DECIMAL ±.001							
FOUR PLACE DECIMAL ±.0002							
INTERPRET GEOMETRIC						KDeb Manufacturing Corporation	
TOLERANCING PER:ANSI/ASME Y14M-1994							
<b>PROPRIETARY AND CONFIDENTIAL</b>							
THE INFORMATION CONTAINED IN THIS DRAWING IS THE SOLE PROPERTY OF BATESVILLE TOOLING & DESIGN, INC. ANY REPRODUCTION IN PART OR AS A WHOLE WITHOUT THE WRITTEN PERMISSION OF BATESVILLE TOOLING & DESIGN, INC IS PROHIBITED.		BTD JOB #				SIZE	
						PART NAME:	
						REV	
						1	
						1 ML KDeb 96 Deep Well Plate	
						SHEET 1 OF 1	
</							