

## YIPF5 Polyclonal Antibody

Catalog number: 12931-1-AP

Size: 67 µg/150 µl

Source: Rabbit

Isotype: IgG

Synonyms:

YIPF5; DKFZp313L2216,

FinGER5, Protein YIPF5, SB140,

SMAP 5, SMAP5,

UNQ3123/PRO10275, Yip1

domain family, member 5, YIP1

family member 5, YIP1A, YIPF5,

YPT interacting protein 1 A



human liver tissue were subjected to SDS PAGE followed by western blot with 12931-1-AP(YIPF5 antibody) at dilution of 1:500

### Applications

Tested applications:	ELISA, WB
Species specificity:	Human,Mouse,Rat; other species not tested.
Calculated YIPF5 MW:	257aa,28kd
Observed YIPF5 MW:	28kd
Positive WB detected in	Human liver tissue, human heart tissue
Recommended dilution:	WB: 1:400 - 1:1600

Application key: WB = Western blotting, IHC = Immunohistochemistry, IF = Immunofluorescence, IP = Immunoprecipitation

### Immunogen information

Immunogen:	Ag3802
GenBank accession number:	BC024737
Gene ID (NCBI):	81555
Full name:	Yip1 domain family, member 5

### Product information

Purification method:	Antigen affinity purification
Storage:	PBS with 0.1% sodium azide and 50% glycerol pH 7.3. Store at -20°C.