

TRIM17 Polyclonal Antibody

Catalog number: 13663-1-AP


Size: 20 µg/150 µl

Source: Rabbit

Isotype: IgG

Synonyms:

TRIM17; RBCC; RNF16; terf



HeLa cells were subjected to SDS PAGE followed by western blot with 13663-1-AP (TRIM17 Antibody) at dilution of 1:600

Applications

Tested applications:	ELISA, WB
Species specificity:	Human, mouse, rat; other species not tested.
Positive WB detected in:	HeLa cells, COLO 320 cells
Calculated TRIM17 MW:	477aa, 54kd
Observed TRIM17 MW:	54kd
Recommended dilution:	WB: 1:200-1:2000

Application key: WB = Western blotting, IHC = Immunohistochemistry, IF = Immunofluorescence, IP = Immunoprecipitation, FC = Flow cytometry

Immunogen information

Immunogen:	Ag4572
GenBank accession number:	BC033788
Gene ID (NCBI):	51127
Full name:	Tripartite motif-containing 17

Product information

Purification method:	Antigen affinity purification
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3. Store at -20°C.