

NUDT2 Polyclonal Antibody

Catalog number: 10484-1-AP


Size: 30 µg/150 µl

Source: Rabbit


Isotype: IgG

Synonyms:

NUDT2; APAH1; MGC10404



mouse liver tissue were subjected to SDS PAGE followed by western blot with 10484-1-AP (NUDT2 Antibody) at dilution of 1:1500



Immunohistochemistry of paraffin-embedded human breast cancer using 10484-1-AP (NUDT2 Antibody) at Dilution 1:100 (under 25x lens)

Applications

Tested applications:	ELISA, WB, IHC
Species specificity:	Human, mouse, rat; other species not tested.
Positive WB detected in:	Mouse liver tissue, HEK-293 cells, mouse spleen tissue, mouse colon tissue
Calculated NUDT2 MW:	17kd
Observed NUDT2 MW:	17kd
Positive IHC detected in:	Human breast cancer
Recommended dilution:	WB: 1:500-1:5000 IHC: 1:50-1:200

Application key: WB = Western blotting, IHC = Immunohistochemistry, IF = Immunofluorescence, IP = Immunoprecipitation, FC = Flow cytometry

Immunogen information

Immunogen:	Ag0743
GenBank accession number:	BC004926
Gene ID (NCBI):	318
Full name:	Nudix (nucleoside diphosphate linked moiety X)-type motif 2

Product information

Purification method:	Antigen affinity purification
Storage:	PBS with 0.1% sodium azide and 50% glycerol pH 7.3. Store at -20°C.