

# Balance Buying Guide

## MS-TS, ML-T, and ME-TE



METTLER TOLEDO recently released new basic and entry level balances: the MS-TS, ML-T, and ME-TE. Each boast a unique set of features and value, and together provide more solutions than ever before. Learn more about the families below, along with the previous generation models as reference to assist you in your selection.

### Need Five Star Performance? ★ ★ ★ ★ ★

METTLER TOLEDO Excellence Balances (not shown here) can support the highest needs for safety, efficiency, quality assurance, and ease of compliance. Please visit [www.thomassci.com](http://www.thomassci.com) or contact your local representative for details on anything you see here as well as our full lineup of Excellence balances and more advanced weighing solutions.


## The New Generation

### MS-TS Balances

#### Sophistication Made Simple

The MS-TS packs high-end features into an easy-to-use instrument. An extra-large 7" color touchscreen display, full-metal housing, and quick-lock draftshield make this a sophisticated balance that's dependable for the everyday user.

- ▶ Analytical
- ▶ Precision




### ML-T Balances

#### Big Features Go Compact

The ML-T balance shares the solid performance of the MS-TS, but with a smaller touchscreen and footprint. Internal calibration, on-screen guidance and a metal bottom housing make it a strong and reliable addition to any lab.

- ▶ Analytical
- ▶ Precision




### ME-TE Balances

#### For Your Essential Needs

The ME-TE balance is an affordable option for those seeking the essentials. A touchscreen and built-in applications combined with a standard draft shield and external calibration make it a quality balance at a friendly price point.

- ▶ Analytical
- ▶ Precision



<b>Performance</b>	★	★	★	★	☆	★	★	★	★	☆	★	★	★	☆	☆
<b>Ease of Use</b>	★	★	★	★	☆	★	★	★	☆	☆	★	★	☆	☆	☆
<b>Durability</b>	★	★	★	★	☆	★	★	★	☆	☆	★	★	☆	☆	☆
<b>Security</b>	★	★	★	★	☆	★	★	★	★	☆	★	★	☆	☆	☆

## The Previous Generation (Phased Out 2015)

### MS Balances

#### Long-Lasting Quality

Looking to replace your METTLER TOLEDO MS Balance? The new MS-TS and ML-T are the closest replacements to this popular product family.



### ML Balances

#### Basic Weighing

Looking to replace your METTLER TOLEDO ML Balance? For ML Analytical balances, the ML-T is your closest replacement. Precision (toploading) owners should check out either the ML-T or the ME-TE (external calibration only).



For pricing and more info, visit [www.thomassci.com](http://www.thomassci.com)

**METTLER TOLEDO**