

SOX6 Polyclonal Antibody

Catalog number: 14010-1-AP


Size: 34 µg/150 µl

Source: Rabbit

Isotype: IgG

Synonyms:

SOX6; SOX6, Transcription
factor SOX 6, SOXD, HSSOX6



Jurkat cells were subjected to SDS
PAGE followed by western blot with
14010-1-AP(SOX6 antibody) at
dilution of 1:500

Applications

Tested applications:	ELISA, WB
Species specificity:	Human, mouse, rat; other species not tested.
Calculated SOX6 MW:	92 kDa
Observed SOX6 MW:	90kd
Positive WB detected in	Jurkat cells
Recommended dilution:	WB: 1:500-1:5000

Application key: WB = Western blotting, IHC = Immunohistochemistry, IF = Immunofluorescence, IP = Immunoprecipitation

Immunogen information

Immunogen:	Ag5111
GenBank accession number:	BC047064
Gene ID (NCBI):	55553
Full name:	SRY (sex determining region Y)-box 6

Product information

Purification method:	Antigen affinity purification
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3. Store at -20°C.