

|               |  |
|---------------|--|
| Product Name: | ABT-HC-AFS-1720  |
| Description:  | <p>General Purpose Auto Defrost Hydrocarbon Laboratory Freezer 17 Cu. Ft. Capacity</p> <p>One solid swing door, right hinged, not reversible</p> <p>Mechanical thermostat</p> <p>Adjustable operating temperature range: -15°C – -20°C</p> <p>Four adjustable shelves and one wire basket</p> <p>Magnetic door gasket for positive seal</p> <p>Heavy duty door handle</p> <p>Leveling legs</p> <p>Probe access port (3/8")</p> <p>White textured cabinet and door</p> <p>High density urethane foam cabinet and door insulation</p> <p>HFC-free refrigerant (R600a), foam insulation and packaging</p> <p>One year parts and labor warranty, plus an additional four year compressor parts warranty</p> <p>Exterior dimensions: 28"W x 29 3/4"D x 68"H</p> <p>Auto defrost</p> <p>115V, 60Hz, 5 Amps, 1/4 HP</p> <p>CSA listed</p> <p>Shipping weight: 248 lbs.</p> <p>DISCLAIMER: Auto defrost freezers incorporate an electric heater installed on the evaporator that assist in the removal of the frost and ice. Although the heater helps speed up the frost and ice removal, it can raise the temperatures within the freezer.</p> |