



Item No.	Qty	Description	Inventory Code	Material
D1	1	250mL SU Storage Bags	ACAX0012	ULDPE
B1	2	1/4"3/8" (6.4*9.6mm) L-Flex Tubing, Length 300mm	DBBX0003	TPE
C1	1	1/4" (6.4mm) Hose barb non-valved in-line coupling body	EAAX0005	PC
C2	1	MPC coupling plug	EAAX0028	PC
C3	1	1/4" (6.4mm) Hose barb non-valved in-line coupling insert	EAAX0017	PC
C4	1	MPC coupling cap	EAAX0025	PC
F1	4	2.5mm cable tie	ADAA0005	Nylon
F2	2	3.6mm cable tie	ADAA0004	Nylon
F3	2	Robert clip (small)	EDBA0026	PA

Customer Approval		LePure Biotech Co., Ltd.				
Company:		Drawing No. HP0250-3AK				
Signature:		Title: SU Storage Bags				
Date:		Designed By:	Date:	Checked By:	Date:	Edition
		Jijia Guo	2023.03.12	Alan	2023.03.12	Rev01
<small>Design Approval the manufacturing of the system(s) based on this design. By signing, you also assume full</small>		<small>We request that you do not disclose this information to any third party involved in this drawing, in order to protect our customer's confidentiality and intellectual property</small>				