

GEMIN7 Polyclonal Antibody

Catalog number: 10625-1-AP

Size: 37 µg/150 µl

Source: Rabbit

Isotype: IgG

Synonyms:

GEMIN7; SIP3

Applications

Tested applications: ELISA
Species specificity: Human, mouse, rat; other species not tested.
Calculated GEMIN7 MW: 15kd
Observed GEMIN7 MW:

Application key: WB = Western blotting, IHC = Immunohistochemistry, IF = Immunofluorescence, IP = Immunoprecipitation, FC = Flow cytometry

Immunogen information

Immunogen: Ag0976
GenBank accession number: BC007793
Gene ID (NCBI): 79760
Full name: Gem (nuclear organelle) associated protein 7

Product information

Purification method: Antigen affinity purification
Storage: PBS with 0.1% sodium azide and 50% glycerol pH 7.3. Store at -20°C, avoid freeze/thaw cycles.