

TXNDC3 Polyclonal Antibody

Catalog number: 13586-1-AP


Size: 38 µg/150 µl

Source: Rabbit


Isotype: IgG

Synonyms:

TXNDC3; CILD6, NME8, Sptrx
2, SPTRX2, TXNDC3



mouse testis tissue were subjected to SDS PAGE followed by western blot with 13586-1-AP(TXNDC3 antibody) at dilution of 1:1000



IP Result of anti-TXNDC3 (IP:13586-1-AP, 4µg; Detection:13586-1-AP 1:1000) with mouse testis tissue lysate 4000µg.

Applications

Tested applications:	ELISA, WB, IP
Species specificity:	Human,Mouse,Rat; other species not tested.
Calculated TXNDC3 MW:	588aa,67 kDa
Observed TXNDC3 MW:	67 kDa
Positive WB detected in	Mouse testis tissue, human testis tissue
Positive IP detected in	Mouse testis tissue
Recommended dilution:	WB: 1:500-1:5000 IP: 1:500-1:5000

Application key: WB = Western blotting, IHC = Immunohistochemistry, IF = Immunofluorescence, IP = Immunoprecipitation

Immunogen information

Immunogen:	Ag4514
GenBank accession number:	BC036816
Gene ID (NCBI):	51314
Full name:	Thioredoxin domain containing 3 (spermatzoa)

Product information

Purification method:	Antigen affinity purification
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3. Store at -20°C.