



## Technical Data

Color	White
Sizes Available	KIT
Packaging	Boxed
Packed	12/Case
Case Dimensions (in)	20.00 x 13.00 x 21.00
Case Weight (lbs)	55.12
Country of Origin	United States
Regulation	--
Kit Size	50 Person
Case Type	Plastic
Construction	--
Certifications	TAA Compliant, ANSI Z308.1-2009
Product Circularity	--

## PIP®

### Contractor First Aid Kit - 50 Person

- Plastic box is gasketed for weatherproof protection of contents
- Convenient molded carry handle includes slots for wall mounting
- Lightweight kit is ideal for use in the office, vehicle or on the job site
- Labels have first aid instructions in English & Spanish
- Meets ANSI Z308.1-2009 requirements
- Can be customized with your logo

## Performance Data

ANSI First Aid Kit Class (2021)	--
ANSI First Aid Kit Class (2015)	--

## Care Instructions

## SKU Information

SKU	COLOR	SIZE	UPC	UPC (INNER)	UPC (CASE)
299-13293	White	KIT	899558001788		20616314847732

**Thomas No. 24A00F983**

PROTECTIVE INDUSTRIAL PRODUCTS, INC. | BRINGING THE BEST OF THE WORLD TO YOU®