

CLEC10A/CD301 Polyclonal Antibody

Catalog number: 13590-1-AP


Size: 27 µg/150 µl

Source: Rabbit

Isotype: IgG

Synonyms:

CLEC10A; CD301, CLEC10A,
CLECSF13, CLECSF14, HML,
HML2, Macrophage lectin 2



BxPC-3 cells were subjected to SDS PAGE followed by western blot with 13590-1-AP (CLEC10A antibody) at dilution of 1:400

Applications

Tested applications:	ELISA, WB
Species specificity:	Human; other species not tested.
Calculated CLEC10A/CD301	316aa, 35 kDa
MW:	
Observed CLEC10A/CD301	32 kDa, 28 kDa
MW:	
Positive WB detected in	BxPC-3 cells
Recommended dilution:	WB: 1:200-1:2000

Application key: WB = Western blotting, IHC = Immunohistochemistry, IF = Immunofluorescence, IP = Immunoprecipitation

Immunogen information

Immunogen:	Ag4490
GenBank accession number:	BC039011
Gene ID (NCBI):	10462
Full name:	C-type lectin domain family 10, member A

Product information

Purification method:	Antigen affinity purification
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3. Store at -20°C.