

CA14 Polyclonal Antibody

Catalog number: 13736-1-AP


Size: 34 µg/150 µl

Source: Rabbit

Isotype: IgG

Synonyms:

CA14;



mouse eye tissue were subjected to SDS PAGE followed by western blot with 13736-1-AP(CA14 Antibody) at dilution of 1:1200

Applications

Tested applications:	ELISA, WB
Species specificity:	Human, mouse, rat; other species not tested.
Positive WB detected in:	Mouse eye tissue
Calculated CA14 MW:	337aa,38kd
Observed CA14 MW:	45-50kd
Recommended dilution:	WB: 1:500-1:5000

Application key: WB = Western blotting, IHC = Immunohistochemistry, IF = Immunofluorescence, IP = Immunoprecipitation, FC = Flow cytometry

Immunogen information

Immunogen:	Ag4673
GenBank accession number:	BC034412
Gene ID (NCBI):	23632
Full name:	Carbonic anhydrase XIV

Product information

Purification method:	Antigen affinity purification
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3. Store at -20°C.