



TheSharpChoice[®]

NMSDC Certified Diversity Supplier

RELI[®] Safety Blood Collection Sets

Reli[®] Safety Blood Collection Sets are OSHA compliant and are available either with or without pre-attached tube holders. All Reli[®] products are manufactured to strict quality standards and contain unique features to improve user convenience, safety, and performance.

Pictured: Reli[®] Safety Blood Collection Set (GSBCS23G-7T)

- Not made with natural rubber latex
- Blood flashback seen clearly through tubing
- Audible click when needle is in final lock position
- Siliconized short bevel needle for smooth insertion
- Flexible, color-coded wing for improved tactile feel and easy gauge identification
- Features single-handed safety shield activation for convenience and ease of use
- Activation of safety shield while needle remains in vein for reduced risk of accidental needle stick injury
- Computerized grinding and 100% visual camera inspection ensures sharp, consistent, quality needles

BRAND (Company)	Color-Coding for Gauge Identification	Audible Click when Shield is Activated	Final Lock Position	Single-Handed In-Vein Needle Retraction
Reli [®] Safety Blood Collection Sets (MYCO Medical)	✓	✓	✓	✓
Vacutainer [®] Safety-Lok [™] Blood Collection Sets (BD)	✓	✓	✓	✓
Vacutainer [®] UltraTouch [™] Push Button Blood Collection Sets (BD)	✓	✓	✓	✓
Vacurette [®] Blood Collection Sets (Greiner Bio-One)	✓	✓	✓	✓
Surshield [®] Safety Blood Collection Set (Terumo)	✓	✓	✓	
K-Shield Advantage [®] Winged Blood Collection Set (Kawasumi)	✓	✓	✓	✓

Reli® Safety Blood Collection Sets (Sterile)

	MYCO CODE	NEEDLE GAUGE	TUBE LENGTH	COLOR	UNIT OF MEASURE
WITHOUT TUBE HOLDER	GSBCS21G-3	21G x 3/4"	10 cm (3.9")	Green	100/bx
	GSBCS21G-7	21G x 3/4"	19 cm (7.4")	Green	100/bx
	GSBCS21G-12	21G x 3/4"	30 cm (11.8")	Green	100/bx
	GSBCS23G-3	23G x 3/4"	10 cm (3.9")	Blue	100/bx
	GSBCS23G-7	23G x 3/4"	19 cm (7.4")	Blue	100/bx
	GSBCS23G-12	23G x 3/4"	30 cm (11.8")	Blue	100/bx
	GSBCS25G-3	25G x 3/4"	10 cm (3.9")	Orange	100/bx
	GSBCS25G-7	25G x 3/4"	19 cm (7.4")	Orange	100/bx
	GSBCS25G-12	25G x 3/4"	30 cm (11.8")	Orange	100/bx
WITH TUBE HOLDER	GSBCS21G-3T	21G x 3/4"	10 cm (3.9")	Green	50/bx
	GSBCS21G-7T	21G x 3/4"	19 cm (7.4")	Green	50/bx
	GSBCS21G-12T	21G x 3/4"	30 cm (11.8")	Green	50/bx
	GSBCS23G-3T	23G x 3/4"	10 cm (3.9")	Blue	50/bx
	GSBCS23G-7T	23G x 3/4"	19 cm (7.4")	Blue	50/bx
	GSBCS23G-12T	23G x 3/4"	30 cm (11.8")	Blue	50/bx
	GSBCS25G-3T	25G x 3/4"	10 cm (3.9")	Orange	50/bx
	GSBCS25G-7T	25G x 3/4"	19 cm (7.4")	Orange	50/bx
	GSBCS25G-12T	25G x 3/4"	30 cm (11.8")	Orange	50/bx

As a leading NMSDC Diversity Certified Supplier, MYCO Medical is committed to “Bringing Reliability, Innovation and Value” to the Healthcare industry with an emphasis on “physician preference” surgical devices and disposables.

To order, learn more or receive complimentary evaluation samples:

Please call MYCO Medical at (800) 454-6926 or (919) 460-2535 or email us at info@mycomedical.com

See the complete line of MYCO innovative disposable and reusable medical devices at www.mycomedical.com



2015 Production Drive
Apex, NC 27539 USA

info@mycomedical.com
www.mycomedical.com

Some medical devices are not licensed for sale in Canada.

© 2019 MYCO Medical, Inc. All rights reserved.
LIT-00547 Rev 0 Eif 01.28.2019