

15900

## **FluidGard® 120 Procedure Mask**

- Complies with ASTM F2100 Level 2 Requirements
- Specially designed inner facing improves comfort
- Integrated layers work together to maximize breathability
- Soft comfortable earloops provide easier donning



| Thomas No. | Mfr. No. | Description  |
|------------|----------|--|
| 22A00C291  | 15900    | Procedure Mask, FluidGard® 120, Blue, ASTM Level 2, Non-Sterile BFE>99.0 |

### **Specifications**

|                   |                               |
|-------------------|-------------------------------|
| Product           | FluidGard® 120 Procedure Mask |
| Color             | Blue                          |
| Fluid Resistance  | 120 mm Hg                     |
| Filtration - BFE  | ≥99.0                         |
| Filtration - PFE  | ≥98.9                         |
| Delta P           | ≤3.3                          |
| Flame Spread      | Class 1                       |
| Attachment Design | Earloop                       |
| Mask Construction | Sonic Welded                  |
| Packaging         | 50/Box, 10/Boxes, 500/Case    |

This product is NOT made with natural rubber latex



### **Complies with ASTM F2100**

|   |                |
|---|----------------|
| <b>Fluid Resistance (mmHg):</b>               | <b>120</b>     |
| <b>BFE (%) @ 3.0 µm:</b>                      | <b>≥ 98%</b>   |
| <b>PFE (%) @ 0.1 µm:</b>                      | <b>≥ 98%</b>   |
| <b>ΔP (mm H<sub>2</sub>O/cm<sup>2</sup>):</b> | <b>&lt; 5</b>  |
| <b>Flame Spread:</b>                          | <b>Class 1</b> |