

TRBV5-4 Polyclonal Antibody

Catalog number: 13524-1-AP

Size: 24 µg/150 µl


Source: Rabbit

Isotype: IgG

Synonyms:

TRBV5-4; TCRBV5S4;

TCRBV5S6A3N2T; TRBV54



U-937 cells were subjected to SDS PAGE followed by western blot with 13524-1-AP (TRBV5-4 Antibody) at dilution of 1:800

Applications

Tested applications:	ELISA, WB
Species specificity:	Human; other species not tested.
Positive WB detected in:	U-937 cells
Calculated TRBV5-4 MW:	308aa, 34kd
Observed TRBV5-4 MW:	30-34kd
Recommended dilution:	WB: 1:200-1:2000

Application key: WB = Western blotting, IHC = Immunohistochemistry, IF = Immunofluorescence, IP = Immunoprecipitation, FC = Flow cytometry

Immunogen information

Immunogen:	Ag4337
GenBank accession number:	BC028083
Gene ID (NCBI):	28611
Full name:	T cell receptor beta variable 5-4

Product information

Purification method:	Antigen affinity purification
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3. Store at -20°C.