

Glass Sprayers

Provides uniformly fine mist or spray, required for paper chromatographic determination

- Conveniently located vent permits instantaneous control of spray
- Borosilicate glass tube has permanently attached, molded phenolic screw cap (size No. 24) with polyvinyl liner which fits commercial narrow mouth bottles
- Complete with spray tube, screw cap and 240 mL clear glass bottle (Thomas No. 1702F61)
- Takes 1/4 inch i.d. rubber tubing

| Thomas No. | Description | Count | Each |
|------------|--------------------|-------|------|
| 2753J10 | Glass Spray Bottle | 6 | 1 |

