




labForce[®] square media bottles are made of high-quality PET material through precision molds and automated equipment, with high transparency and high durability. These bottles are ideal for long-term storage of aqueous solutions, such as the storage of serum, culture medium, pancreatin, reagents and other products.



PRODUCT FEATURES:

- *Bottle material: PET, cap material: HDPE
- *Bottle with precise scale
- *The bottle is square, easy to hold and convenient to store and transport.
- *Available in 30ml, 60ml, 125ml, 250ml, 500ml, 1000ml
- *Rigorous leak testing to ensure no leaks
- *RNase\DNase free, pyrogen free, endotoxin free
- *Can be stored frozen at temperatures as low as -135° C
- *Sterilized by E-beam irradiation
- *Packed in shrink-wrapped tray modules.

Thomas No.	Description	Qty/List
21A00Q104	labForce [®] Sterile Square PET Media Bottles 30ml (40/PK, 7 PKS/CS)	280/Cs \$350.10
21A00Q105	labForce [®] Sterile Square PET Media Bottles 60ml (40/PK, 5 PKS/CS)	200/Cs \$357.70
21A00Q099	labForce [®] Sterile Square PET Media Bottles 125ml (24/PK, 4 PKS/CS)	96/Cs \$271.20
21A00Q100	labForce [®] Sterile Square PET Media Bottles 250ml (30/PK, 2 PKS/CS)	60/Cs \$259.50
21A00Q101	labForce [®] Sterile Square PET Media Bottles 500ml (20/PK, 2 PKS/CS)	40/Cs \$225.45
21A00Q102	labForce [®] Sterile Square PET Media Bottles 1000ml (12/PK, 2 PKS/CS)	24/Cs \$129.90