

PNMA1 Polyclonal Antibody

Catalog number: 13631-1-AP


Size: 37 µg/150 µl

Source: Rabbit


Isotype: IgG

Synonyms:

PNMA1; MA1



mouse testis tissue were subjected to SDS PAGE followed by western blot with 13631-1-AP(PNMA1 Antibody) at dilution of 1:1000



Immunohistochemistry of paraffin-embedded human gliomas using 13631-1-AP(PNMA1 Antibody) at Dilution 1:100 (under 40x lens)

Applications

Tested applications:	ELISA, WB, IHC
Species specificity:	Human, mouse, rat; other species not tested.
Positive WB detected in:	Mouse testis tissue
Calculated PNMA1 MW:	353aa,40kd
Observed PNMA1 MW:	37kd
Positive IHC detected in:	Human gliomas
Recommended dilution:	WB: 1:500-1:5000 IHC: 1:50-1:200

Application key: WB = Western blotting, IHC = Immunohistochemistry, IF = Immunofluorescence, IP = Immunoprecipitation, FC = Flow cytometry

Immunogen information

Immunogen:	Ag4549
GenBank accession number:	BC039577
Gene ID (NCBI):	9240
Full name:	Paraneoplastic antigen MA1

Product information

Purification method:	Antigen affinity purification
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3. Store at -20°C.