

GEB® 1000µL Universal Filter Tips, Extended, Clear, Sterile, Rack Pack, 96 Tips/Rack, 10 Racks/Pack, 5 Packs/Case (Product # FTL1000-R-NS)

GEB® Filter Barrier tips are ideal for DNA amplification and microbiological applications where trace contamination can completely invalidate results. The filter prevents cross-contamination from sample to sample by creating a physical barrier in the tip to inhibit and block aerosols. Every polyethylene filter, free from DNA inhibitors, is inserted into the tip robotically to ensure freedom from nuclease contamination. All GEB filter tips comply with USP 3127 (Heavy Metal Contamination) and are certified to be RNase, DNase, and endotoxin free. They are packaged in sealed-bottom racks with lids to prevent contamination and are available in sterilized form.

Certified to be RNase-/DNase-free


Nonpyrogenic

Sterile

Product details at-a-glance

<b>Qty./Pack.:</b>	960 pieces tips
<b>Qty./Case.:</b>	4800 pieces tips
<b>Color:</b>	Clear
<b>Sterile:</b>	Yes
<b>Packaging Format:</b>	96 Tips/Rack, 10 Racks/Pack, 5 Packs/Case
<b>Case dimension:</b>	495*285*650mm

Drawing:



Prepared by: Anne June 16, 2020 Approve: WeiMing Yu June 16, 2020