

## PTGR2 Polyclonal Antibody

Catalog number: 14164-1-AP

Size: 56 µg/150 µl

Source: Rabbit


Isotype: IgG

Synonyms:

PTGR2; DKFZp686P10120;

FLJ39091; FLJ99229; PGR2;

ZADH1



mouse kidney tissue were subjected to SDS PAGE followed by western blot with 14164-1-AP(PTGR2 Antibody) at dilution of 1:500

### Applications

Tested applications:	ELISA, WB
Species specificity:	Human, mouse, rat; other species not tested.
Positive WB detected in:	Mouse kidney tissue
Calculated PTGR2 MW:	351aa,39kd
Observed PTGR2 MW:	38.5kd
Recommended dilution:	WB: 1:200-1:2000

Application key: WB = Western blotting, IHC = Immunohistochemistry, IF = Immunofluorescence, IP =

Immunoprecipitation, FC = Flow cytometry

### Immunogen information

Immunogen:	Ag5352
GenBank accession number:	BC059364
Gene ID (NCBI):	145482
Full name:	Prostaglandin reductase 2

### Product information

Purification method:	Antigen affinity purification
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3. Store at -20°C.