


Non-Invasive Temperature Readings for the Workforce

Peel and Stick Thermometers Offer Non-Invasive Method for Safety Protocol



H-B DURAC® Liquid Crystal Disposable Body Thermometer

Fast and Reliable Method for Monitoring Body Temperature

- Non-invasive body temperature measurement ideal for workforce personnel
- Color-coded readings allow for quick evaluation
- Can be self-applied to maintain social distancing during temperature testing
- Easy to use - simply peel off backing and stick to forehead or hand
- Latex-free
- Utilizes Thermocratic Liquid Crystal technology
- Disposable
- Instructions included



Catalog No.	Range	Divisions	Accuracy	Length	Quantity/Pack
B62010-0000	94.6/104.6°F	2	±1°F	45mm (1.8")	100

