




## Color Wheel

Discover the principles of light and color with this dynamic color wheel. The 8-inch disc features vibrant color segments that blend into white when spun, highlighting the science behind color and additive color mixing. This demonstrates that white light contains the complete visible spectrum through the combination of individual color wavelengths.

The apparatus is made of ABS plastic and features a smooth belt-drive mechanism and ergonomic hand crank for reliable classroom use. Perfect for physics and optics curricula, this tool brings color theory concepts to life through hands-on experimentation. Includes access to a downloadable activity guide.



Hand-crank operation for interactive learning.

- + Demonstrates additive color theory principles
- + Belt-driven system ensures smooth rotation
- + Ideal for classrooms, science museums and curious minds
- + Downloadable activity guide

### Color Wheel

| Item #  | Description | Quantity Per Pack |
|---------|-------------|-------------------|
| COLRWHL | Color Wheel | 1                 |



**WARNING:** Can expose you to acrylonitrile and styrene, carcinogens, and 1,3-butadiene, a carcinogen and reproductive toxicant. See [www.P65Warnings.ca.gov](http://www.P65Warnings.ca.gov).



1001 Technology Way | Libertyville, IL 60048

[www.unitedsci.com](http://www.unitedsci.com) | 847-336-7556



**United  
Scientific®**