

Recombinant Influenza Hemagglutinin (H5N1) (A/Vietnam/1203/2004)

Catalog #: 1501-001

Lot #: 1509005


Description: Recombinant, Influenza Hemagglutinin Full-Length Glycoprotein from virus strain A/Vietnam/1203/2004 (H5N1). Recombinant HA (rHA) is supplied as an affinity purified protein. rHA is produced in Sf9 insect cells using baculovirus for expression and is purified by FPLC.

Storage: 2-3 weeks at -20°C, -80°C long term

Size: 50 µg of protein is supplied at a concentration of 0.362 mg/mL (by BCA) in Tris buffered saline and 0.01% non-ionic surfactant. The theoretical molecular weight of the protein is ~64 kDa, without glycosylation. Because of the highly glycosylated nature of this protein, migration in an SDS-PAGE gel is slowed resulting in broad, diffuse bands representing differing glycosylation forms.

Relevance: Recombinant hemagglutinin glycoprotein provides a control protein for immunoassays and a tool to enhance Orthomyxovirus research.


SDS-PAGE and Western Blot Detection



(A) SDS-PAGE and stain demonstrating 1 µg and 5 µg (lanes 1 and 2, respectively) of rHA protein under denaturing and reducing conditions. (B) Western blot detection of rHA at 200 ng when detected with anti-H5 chicken antiserum (1:1,000 dilution) and an anti-chicken AP conjugate followed by substrate.

Hemagglutination using Turkey Red Blood Cells:

rHA lot 1509005 showed HA titer of 1:8192.



For additional influenza products, please visit:
<http://ibtbioservices.com/index.php/product-and-reagents/influenza-products>

Intended for research use only. Not for human, therapeutic, or diagnostic applications.
 The buyer cannot sell or otherwise transfer this product for Commercial Purposes without written approval of Integrated BioTherapeutics, Inc.

Copyright 2015. Integrated BioTherapeutics, Inc. All rights reserved