

## RARRES2 Polyclonal Antibody

Catalog number: 10216-1-AP

Size: 21 µg/150 µl

Source: Rabbit

Isotype: IgG

Synonyms:

RARRES2; RARRES2, Chemerin,

HP10433, RAR responsive

protein TIG2, TIG2

### Applications

Tested applications: ELISA

Species specificity: Human; other species not tested.

Calculated RARRES2 MW: 163aa,19 kDa

Observed RARRES2 MW:

Application key: WB = Western blotting, IHC = Immunohistochemistry, IF = Immunofluorescence, IP = Immunoprecipitation

### Immunogen information

Immunogen: Ag0261

GenBank accession number: BC000069

Gene ID (NCBI): 5919

Full name: Retinoic acid receptor responder (tazarotene induced)  
2

### Product information

Purification method: Antigen affinity purification

Storage: PBS with 0.1% sodium azide and 50% glycerol pH  
7.3. Store at -20°C.