

Efficient PCR Sampling & Storage

Designed for pre- and post-PCR sampling, streamline your PCR workflow with the 32-Place PCR Rack. This versatile rack accommodates 0.2 mL tubes both individually and in strips of 8, as well as 32-place assay plates (8x4)

Space Saving Design

- Small and compact design is stackable, saving valuable space in the freezer
- Right sized when processing smaller batches
- Lid hinge design allows lid to sit at 180 degrees, clear of sample processing, or can be removed entirely

Reusable Convenience

- Transparent lid protects samples and allows for quick assessment of contents
- Embossed numbers and letters on the wells provide easy identification and tracking of samples
- Light-weighted construction made durable for reusability



Hinged Lid
stays clear of sample processing



Embossed Letters & Numbering
on the wells allow for reliable tracking of samples



Made from 100% recyclable Polypropylene

PK10 / CS 20

| Item No. | Tubes mL | No. of Wells | L x W x H in | L x W x H cm |
|----------|-------------|--------------|-----------------|------------------|
| 120889 | 0.2 | 32 | 3.8 x 6.9 x 1.2 | 9.8 x 17.5 x 3.1 |

Stores down to -80°C