

PLEKHA3 Polyclonal Antibody

Catalog number: 11292-1-AP

Size: 25 µg/150 µl


Source: Rabbit

Isotype: IgG

Synonyms:

PLEKHA3; FAPP 1, FAPP1,

FLJ20067, PLEKHA3



COLO 320 cells were subjected to SDS PAGE followed by western blot with 11292-1-AP (PLEKHA3 antibody) at dilution of 1:200

Applications

Tested applications:	ELISA, WB
Species specificity:	Human, Mouse, Rat; other species not tested.
Calculated PLEKHA3 MW:	34 kDa
Observed PLEKHA3 MW:	34 kDa
Positive WB detected in	COLO 320 cells, Jurkat cells, K-562 cells, mouse brain tissue
Recommended dilution:	WB: 1:200-1:2000

Application key: WB = Western blotting, IHC = Immunohistochemistry, IF = Immunofluorescence, IP = Immunoprecipitation

Immunogen information

Immunogen:	Ag1851
GenBank accession number:	BC044567
Gene ID (NCBI):	65977
Full name:	Pleckstrin homology domain containing, family A (phosphoinositide binding specific) member 3

Product information

Purification method:	Antigen affinity purification
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3. Store at -20°C.