




UV Flashlight 365 nm

This compact UV flashlight delivers bright 365 nm illumination for convenient on-the-go inspection and fluorescence work.

This UV-A flashlight is ideal for forensics applications or for highlighting fluorescent and blacklight-reactive substances in classroom, laboratory, or field environments. The long-life LED emits at a wavelength of 365 nm (UV-A) and produces a broad beam that illuminates a large area for easy visualization. An attached lanyard provides a secure grip to help prevent drops during use. The built-in rechargeable battery connects to any standard USB port for charging on the go. A USB 2.0 Micro B charging cord is included. For best performance, do not operate the flashlight for more than 8 continuous minutes. Flashlight measures 3 1/4 in long x 3/4 in diameter.

Includes UV-A LED light flashlight with plastic housing and attached lanyard, a rechargeable battery, and a USB 2.0 Micro B charging cord.



UV light demonstration.

- + Bright 365 nm UV-A LED with long service life
- + Rechargeable design with included USB 2.0 Micro B charging cord
- + Compact housing with lanyard and wide-beam output to illuminate a large area

UV flashlight 365 nm

Item #	Description	Quantity Per Pack
UWLIGHT365	UV Flashlight 365 nm	1



WARNING: Can expose you to Lead, a carcinogen and reproductive harm. See www.P65Warnings.ca.gov.



1001 Technology Way | Libertyville, IL 60048

www.unitedsci.com | 847-336-7556



**United
Scientific®**