

## Goat-IgG Polyclonal Antibody

Catalog number: 10285-1-AP

Size: 60 µg/150 µl

Source: Rabbit

Isotype: IgG

Synonyms:

goat IgG

### Applications

**Tested applications:** ELISA  
**Species specificity:** Goat; other species not tested.  
**Calculated Goat-IgG MW:** 155 kDa  
**Observed Goat-IgG MW:**

Application key: WB = Western blotting, IHC = Immunohistochemistry, IF = Immunofluorescence, IP = Immunoprecipitation

### Immunogen information

**Immunogen:** Recombinant Protein  
**GenBank accession number:**  
**Gene ID (NCBI):**  
**Full name:** Undefined

### Product information

**Purification method:** Antigen affinity purification  
**Storage:** PBS with 0.1% sodium azide and 50% glycerol pH 7.3. Store at -20°C.