

HLA-DPB1 Polyclonal Antibody

Catalog number: 11298-1-AP

Size: 20 µg/150 µl

Source: Rabbit


Isotype: IgG

Synonyms:

HLA-DPB1; DPB1, HLA DP1B,

HLA DPB1, HLA-DPB1, MHC

class II antigen DPB1



Immunohistochemistry of paraffin-embedded human spleen tissue slide using 11298-1-AP (HLA-DPB1 Antibody) at dilution of 1:50 (under 40x lens)

Applications

Tested applications:	ELISA, IHC
Species specificity:	Human; other species not tested.
Calculated HLA-DPB1 MW:	29 kDa
Observed HLA-DPB1 MW:	
Positive IHC detected in	Human spleen tissue, human tonsillitis tissue
Recommended dilution:	IHC: 1:20-1:200

Application key: WB = Western blotting, IHC = Immunohistochemistry, IF = Immunofluorescence, IP = Immunoprecipitation

Immunogen information

Immunogen:	Ag1815
GenBank accession number:	BC015000
Gene ID (NCBI):	3115
Full name:	Major histocompatibility complex, class II, DP beta 1

Product information

Purification method:	Antigen affinity purification
Storage:	PBS with 0.1% sodium azide and 50% glycerol pH 7.3. Store at -20°C.