



Recommended Uses



General Description

Material	Synthetic Vinyl	Type	Non-Sterile
Color	Blue	Usage	Single Use
Shape	Ambidextrous	Interior	Powder Free
Sizes	S - XL	Exterior	Smooth
Cuff	Beaded	Process	Polymer Coated

Product Codes

Small	20A00U985
Medium	20A00V254
Large	20A00U987
X-Large	20A00U988

Packaging

Box	100 Gloves Per Box
Box Dims	230 x 120 x 60 mm
Case	10 Boxes Per Case
Case Dims	315 x 245 x 250 mm

Physical Properties

Tensile Strength	14 Mpa (min)
Tensile Strength (Aged)	14 Mpa (min)
Powder Content	< 2 mg
FDA Listing	510(k) Status
Elongation	350%
Elongation (Aged)	350%
Protein Content	N/A
AQL Level	2.5

Dimensions

Length	235 ± 5 mm		
Width		Weight	Thickness
Small	85 ± 5 mm	Small	4.5 ± 0.3 gm
Medium	95 ± 5 mm	Medium	5.0 ± 0.3 gm
Large	105 ± 5 mm	Large	5.3 ± 0.3 gm
X-Large	115 ± 5 mm	X-Large	5.6 ± 0.3 gm
		Cuff	3 ± 1 mil
		Palm	3 ± 1 mil
		Finger	3 ± 1 mil

Barcodes

 6 97383 97413 6 S Box	 6 97383 97414 3 M Box	 6 97383 97415 0 L Box	 6 97383 97416 7 XL Box
 6 97383 97418 1 S Case	 6 97383 97419 8 M Case	 6 97383 97420 4 L Case	 6 97383 97421 1 XL Case