

MAK Polyclonal Antibody

Catalog number: 13638-1-AP


Size: 29 µg/150 µl

Source: Rabbit

Isotype: IgG

Synonyms:

MAK; dJ417M14.2



PC-3 cells were subjected to SDS PAGE followed by western blot with 13638-1-AP(MAK Antibody) at dilution of 1:600

Applications

Tested applications:	ELISA, WB
Species specificity:	Human, mouse, rat; other species not tested.
Positive WB detected in:	PC-3 cells, HeLa cells
Calculated MAK MW:	623aa,71kd
Observed MAK MW:	60-62kd
Recommended dilution:	WB: 1:200-1:2000

Application key: WB = Western blotting, IHC = Immunohistochemistry, IF = Immunofluorescence, IP = Immunoprecipitation, FC = Flow cytometry

Immunogen information

Immunogen:	Ag4571
GenBank accession number:	BC039825
Gene ID (NCBI):	4117
Full name:	Male germ cell-associated kinase

Product information

Purification method:	Antigen affinity purification
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3. Store at -20°C.